




Beamline Staff:

Dr. Richard Gillilan
Dr. Qingqui Huang
Dr. Durgesh Rai

Contributors:

Dr. Nozomi Ando
Robert Miller

High-Pressure Applications:

'Anything a molecule does
changes it's volume'

- Deep Life
- Food Processing
- Pharmaceuticals
- General Biophysical Tool

Beamline Details:

- New Eiger4M Detector
- Highest flux: 7×10^{12} ph/s
- q -range ($\text{\AA}^{-1}$): 0.005-0.7
- 0-100 MPa pressure range
- Real time data readout

